



GranToxiLux®

(Updated 04/09/09)

Overview

GranToxiLux® is OncoImmunitin, Inc.'s single cell-based fluorogenic cytotoxicity assay for the measurement of granzyme B activity in live target cells following the successful delivery of a lethal hit by cytotoxic lymphocytes. Granzyme B exists in lysosomal granules in an inactive form in effector, e.g., CTL and NK, cells. When degranulation of effectors is induced in the presence of target cells, Granzyme B as well as other granule contents such as perforin are taken up by the latter with Granzyme B believed to be the first active enzyme delivered to target cells. Thus, measurement of granzyme B activity inside target cells provides an extremely early quantitative assessment of cell-mediated cellular cytotoxicity.

GranToxiLux® is similar to CyToxiLux® PLUS! and PanToxiLux™, OncoImmunitin, Inc.'s other single cell cytotoxicity assay kits, with the difference being the cell permeable, fluorogenic substrate: CyToxiLux® PLUS! is designed to detect downstream caspase activity only and PanToxiLux™ detects both granzyme B and upstream caspase activities. As with the other two kits, GranToxiLux® can be used for selection of antibodies operating via an antibody-dependent cellular cytotoxicity (ADCC) mechanism in both low and high throughput screening (HTS) modes.

Advantages of GranToxiLux®, CyToxiLux®, and PanToxiLux™ over other cytotoxicity assays, e.g., ⁵¹Cr release, LDH release, and PI, include: (1) cytotoxicity is measured as a fundamental biochemical pathway leading to cell death (cleavage of a cell permeable fluorogenic substrate) rather than merely as the loss of plasma membrane permeability and its sequelae, (2) sensitivity is enhanced such that relatively weak CTL responses against subdominant epitopes are detectable (3) rapidity (Effector:Target coincubation times between 0.3 and 2 hours), (4) measurement of cell death can be carried out exclusively in target cell populations by flow cytometry or fluorescence microscopy, and (5) when combined with immunophenotypic analyses and multiparameter flow cytometry, cytotoxic lymphocyte-mediated killing of primary host target cells as well as the physiology and fate of effector cells can be directly visualized and monitored.

Target cells are fluorescently labeled (red) and then coincubated with cytotoxic effector cells in the presence of a fluorogenic granzyme B substrate. Following incubation and washing, samples may be analyzed by flow cytometry. Cleavage of the substrate results in increased green fluorescence in dying cells. Real-time imaging can also be carried out with confocal microscopy.

Please read this entire protocol before commencing assay!

Components supplied in GranToxiLux® kit (sufficient for 50 assays)

Vial GS (3 vials) = Granzyme B Substrate solution

Vial TFL2 (1 vial) = Target cell marker for use with single laser instruments

Vial TFL4 (1 vial) = Target cell marker for use with dual laser instruments

TFL dilution medium (2 small Eppendorf tubes) = Resuspension medium

for Vials TFL2 and TFL4

Wash Buffer bottle (1 bottle)

Components supplied by user

Effector cells

Target cells

In order to eliminate target cells that have died prior to addition of effector cells, e.g., due to a freeze/thaw cycle, please contact OncoImmunitin, Inc. Please note: labeling of cell surfaces with fluorescent probes is **not** recommended as this may interfere with effector-target interactions.

Format: The assay may be performed using either 96-well plates or polypropylene microcentrifuge tubes.

Reconstitution of TFL2 or TFL4: Add 25 µl from one of the TFL dilution medium tubes to **Vial TFL2** or **Vial TFL4**. (Once reconstituted, **Vial TFL2** or **TFL4** should be stored at -20°C.) **TFL2** is for use with flow cytometers equipped with only Ar ion (488 nm) lasers; measurements require compensation (please see **FLOW CYTOMETRY** instructions below). **TFL4** is for use with flow cytometers equipped with both Ar ion (488 nm) and red He-Ne (633 nm) lasers; measurements do not require compensation. (Neither **Target cell** marker is CTO™, the probe used in the **Nature Med.** reference below; rather, they are dyes with fluorescence properties more complementary to those of the **Granzyme B Substrate**.)

Medium T = Medium for labeling target cells. 1 µl from either reconstituted **TFL2** or **TFL4** is added to the medium in which target cells had been grown or to a physiologic buffer such as phosphate buffered saline (PBS). Please note: most cells load more efficiently in PBS or in a medium free of serum. If 10% serum is included, the recommended **TFL** dilution is 1:1000 whereas for serum-free buffers, the **Target cell** marker is typically used at 1:3000; however, further dilution may be superior. *Optimal TFL2 or TFL4 concentrations should be determined for individual target cell types as 1:1000 for serum-containing and 1:3000 for serum-free media are merely suggested starting points.*

Washing is defined as addition of the indicated volume of medium/buffer followed by centrifugation and then careful removal of all liquid from tubes or flicking followed by light tamping of plates. Resuspension of pellets should be carried out with gentle pipetting of and/or tapping of tubes with finger. **DO NOT VORTEX.**

Protocol

Preparation of Target cells

1. Suspend Target cells (suspension or trypsinized adherent cells) in **Medium T** at 2×10^6 cells/ml. (This is a suggested concentration. Lower numbers can and are routinely used. The critical point is to be able to collect 5,000-10,000 Target cells for analysis.) If targets are to be pulsed with sensitizers, e.g., peptides, the latter should be added at this stage. (If peptide solubility requires use of an organic cosolvent, the latter should not exceed 0.3% (v/v) and a vehicle control tube/well should be included.) Incubate at 37°C for 0.25-1.0 hour. Optimal time should be determined for individual cell types (and sensitizers). Please note: for loading of PBLs, 1:1000 dilution of **TFL2** or **TFL4** in PBS for 15' at 37 °C is recommended. If sensitizers are to be used with PBLs, the recommendation is to add **TFL2** or **TFL4** for the final 15' of sensitizer exposure.
2. During this time, prepare Effector cells (see below).
3. Wash **Target cells** 2 times with at least a 10-fold excess volume of physiologic buffer/medium per wash.
4. Resuspend labeled **Target cells** at 1×10^6 cells/ml in **Wash Buffer**. (Depending on the experimental design, lower numbers of **Target cells** may be used.)
5. Dispense 100 μ l of **Target cell** suspension to each assay well or tube.

Preparation of Effector cells

1. Prepare Effector cells at the appropriate concentration in **Wash Buffer**. For example, for a final Effector to **Target** ratio of 5:1, suspend Effector cells at 5×10^6 cells/ml.

Coincubation of Target and Effector cells

1. Add 100 μ l of Effector cell suspension to each well (or tube) containing **Target cells** except at least two wells, and add 100 μ l of Effector cell suspension to at least two wells which do not contain **Target cells**.
2. Add 100 μ l of **Wash Buffer** to wells containing only **Targets** and only Effectors to bring all samples to a final volume of 200 μ l.
3. Centrifuge all samples, carefully remove medium, and resuspend cell pellets in either 75 μ l **Substrate** from **Vial GS** or **Wash Buffer** for controls (absence of **Substrate** (for Tubes A and C below)). **For ADCC antibodies should be added at this time.**
4. Immediately after resuspension, pellet cells by brief centrifugation (Once the centrifuge reaches the speed normally used for cells, hold for 30 seconds and then stop the instrument.).
5. Incubate at 37°C for the desired time points. Since this assay detects dying cells rather than cells with irreversibly damaged plasma membranes, incubation times for a given cell system should be significantly shorter than with other methodologies. For single point assays, the suggested coincubation time is 1 hour (enduser times range between 30 minutes and 2 hours). Substantially longer times are not recommended.
6. Wash each sample with 200 μ l **Wash Buffer**.
7. Resuspend each sample in **Wash Buffer**, transfer to flow cytometry tubes or leave in plates if a plate reader is to be used, and analyze by flow cytometry.

Summary of samples:

- A. **Target cells** (=Target cells loaded with **TFL2** or **TFL4**)
- B. **Target cells** + **Substrate** from **Vial GS**
- C. Effector cells
- D. Effector cells + **Substrate** from **Vial GS**
- E. **Target cells** + Effector cells + **Substrate** from **Vial GS** (multiple samples)

Flow Cytometry: the following channels are recommended: **APC** for **TFL4**, **PE** for **TFL2** and **FITC** for **CG**. For instruments with other filter sets, please use settings consistent with the following excitation and emission peaks:

TFL4 λ_{ex} : 633 nm, λ_{em} : 657 nm

TFL2 λ_{ex} : 488 nm, λ_{em} : 580 nm

GS λ_{ex} : 488 nm, λ_{em} : 520 nm

For flow cytometers equipped with both Ar ion (488 nm) and red He-Ne (633 nm) lasers (**TFL4** is used as target cell dye):

1. Use samples **A** and **D** to initially set the **APC** and **FITC** channels, resp.: place the peak for cells from sample **A** near 10^3 in the **APC** channel and the peak from sample **D** at ca. 10^1 in the **FITC** channel. Cells in sample B should then be at ca. 10^3 in the **APC** and 10^1 in the **FITC** channels. (Note: Healthy (>95% viable by Trypan Blue) **TFL4**-labeled **Target cells** should appear as a single population. If not, first try decreasing labeling time to 15 minutes. If more than one population still exists, decrease **TFL4** concentration.)
2. Run remaining samples. 5,000-10,000 **Target cells** per sample should be collected for analysis.

For flow cytometers equipped with only an Ar ion (488 nm) laser (**TFL2** is used as target cell dye):

1. Use samples A and D to initially set the **APC** and **FITC** channels, resp. Place the peak from sample **A** near 10^3 in the **APC** channel and the peak from sample **D** at ca. 10^1 in the **FITC** channel. (Note: Healthy (>95% viable by Trypan Blue) **TFL2**-labeled **Target cells** should appear as a single population. If not, first try decreasing labeling time to 15 minutes. If more than one population still exists, decrease **TFL4** concentration.)
2. Use sample **B** to set up **PE** compensation. Dead/dying Effector cells may show a high **FITC** x **PE** population on most single-laser flow cytometers. Compensate **PE** by **FITC** until this population is on the same horizontal axis as viable Effector cells (low **PE**).
3. Run remaining samples. 5,000-10,000 **Target cells** per sample should be collected for analysis.
(Note: channel amplification voltage settings for **FITC** and **PE** after these 3 steps should be similar in magnitude for ease in compensation. Should the **PE** voltage be greater than **FITC** by more than 100, then increase **TFL2** labeling concentration or remove or reduce serum content.)

References for **GranToxiLux®**

1. B.Z. Packard, W.G. Telford, A. Komoriya, and P.A. Henkart. Granzyme B activity in target cells detects attack by cytotoxic lymphocytes. **J. Immunol.** *179*:3812-3820 (2007).
2. A. Kinter, J. McNally, L. Riggins, R. Jackson, G. Roby, and A.S. Fauci. Suppression of HIV-specific T cell activity by lymph node CD25⁺ regulatory T cells from HIV-infected individuals. **PNAS** *104*:3390-3395 (2007).
3. T. Tarasenko, H.K. Kole, A.W. Chi, M.M. Mentink-Kane, T.A. Wynn, and S. Bolland. T cell-specific deletion of the inositol phosphatase SHIP reveals its role in regulating Th1/Th2 and cytotoxic responses. **PNAS** *104*(27):11382-11387 (2007).
4. M. Carlsten, N.K. Björkström, H. Norell, *et al.* DNAX Accessory Molecule-1 mediated recognition of freshly isolated ovarian carcinoma by resting Natural Killer cells. **Cancer Research** *67*: 1317-1325 (2007).
5. A.L. Kinter, R. Horak, M. Sion, *et al.* CD25⁺ regulatory T Cells isolated from HIV-infected individuals suppress the cytolytic and nonlytic antiviral activity of HIV-specific CD8⁺ T Cells *in vitro*. **AIDS Research and Human Retroviruses** *23*:438-450 (2007).
6. R. Chiarle, C. Martinengo, C. Mastini, C. Ambrogio, V. D'Escamard, G. Forni, and G. Inghirami. The anaplastic lymphoma kinase is an effective oncoantigen for lymphoma vaccination. **Nat. Med.** *14*:676-80 (2008).
7. A. Singh, H. Nie, B. Ghosn, H. Qin, L.H. Kwak, and K. Roy. Efficient modulation of T-cell response by dual-mode, single-carrier delivery of cytokine-targeted siRNA and DNA vaccine to antigen-presenting cells. **Mol. Ther.** *16*:2011-2021 (2008).
8. S.A. Migueles, C.M Osborne, C. Royce *et al.* Lytic granule loading of CD8⁺ T cells is required for HIV-infected cell elimination associated with immune control. **Immunity** *29*:1009-21.

Sample Flow Cytometry Data

Target cells (Jurkat, K562, or Daudi) were incubated in 75 μ l **Granzyme B Substrate** with or without Effector cells (NK or CTL) at a 5:1 Effector:Target ratio for 60 minutes. (CTL assays were carried out using a redirected cytotoxicity protocol. Please contact OncoImmunit, Inc. for details.) Quadrants **R1** (upper left of each panel) represent viable target cells while quadrant **R2** (upper right) represent dying, **Granzyme B**-positive target cells. Effector cells occupy the lower 2 quadrants in Effector + Target samples. The inset **% values** are calculated as $R2/(R1+R2)$.

