



NFL1

(Updated 05/06/09)

Overview

NFL1 is OncoImmunitin, Inc.'s probe for eliminating target cells that have died prior to the start of experiments in which CyToxiLux[®], GranToxiLux[®], or PanToxiLux[™] are used to detect cell-mediated cytotoxicity. Target cells are suspended in a solution containing NFL1 and TFL4, washed and then exposed to a cell permeable substrate with or without effectors. Cells that are positive for NFL1 are gated out; thus, the background of target cells due to low viability resulting from freeze/thaw cycles or suboptimal culturing is eliminated. Please note: both NFL1 and TFL4 are intracellular stains and, therefore, do not interfere with recognition between effectors and targets.

Please read this entire protocol before commencing assay!

NFL1 is supplied as a 1000X solution (for most applications) and should be stored at 4 °C.

Reconstitution of TFL4 (from CyToxiLux[®], GranToxiLux[®], or PanToxiLux[™] kits): Add 25 µl from the TFL4 dilution medium tube to **Vial TFL4**. (Once reconstituted, TFL4 should be stored at -20°C.)

Medium T = Medium for labeling target cells. 1 µl from reconstituted TFL4 and 1 µl from the supplied NFL1 solution are added to the medium in which target cells had been grown or to a physiologic buffer such as phosphate buffered saline (PBS). Please note: most cells load more efficiently in PBS or in a medium free of serum. If 10% serum is included, the recommended TFL4 and NFL1 dilutions are both 1:1000 whereas for serum-free buffers, TFL4 is typically used at 1:3000 but NFL1 dilution may remain at 1:1000; however, further dilution, particularly of the TFL4 may be superior. *Optimal TFL4 concentrations should be determined for individual target cell types as 1:1000 for serum-containing and 1:3000 for serum-free media are merely suggested starting points.* A 1:1000 dilution of NFL1 can be used with most media.

Washing is defined as addition of the indicated volume of medium/buffer followed by centrifugation and then careful removal of all liquid from tubes or flicking followed by light tamping of plates. Resuspension of pellets should be carried out with gentle pipetting of and/or tapping of tubes with finger. **DO NOT VORTEX.**

Protocol

Preparation of Target cells

1. Suspend Target cells (suspension or trypsinized adherent cells) in **Medium T** at 2×10^6 cells/ml and label as **sample A**. (This is a suggested concentration. Lower numbers can and are routinely used. The critical point is to be able to collect 5,000-10,000 Target cells for analysis.) A second tube, labeled **sample Ac**, of approximately 1×10^6 cells in 500 µl of Medium T minus NFL1 should also be incubated under the same conditions. If targets are to be pulsed with sensitizers, e.g., peptides, the latter should be added at this stage. (If peptide solubility requires use of an organic cosolvent, the latter should not exceed 0.3% (v/v) and a vehicle control tube/well should be included.) Incubate at 37°C for 0.25-1.0 hour. Optimal time should be determined for individual cell types (and sensitizers). Please note: for loading of PBLs, 1:1000 dilution of both TFL4 and NFL1 in PBS for 15' at 37 °C is recommended. If sensitizers are to be used with PBLs, the recommendation is to add TFL4 and NFL1 for the final 15' of sensitizer exposure.
2. During this time, prepare Effector cells (see below).
3. Wash **Target cells** 2 times with at least a 10-fold excess volume of physiologic buffer/medium per wash.
4. Resuspend labeled **Target cells** at 1×10^6 cells/ml in **Wash Buffer**. (Depending on the experimental design, lower numbers of **Target cells** may used.) Set aside **sample Ac** until flow cytometry measurements commence and continue setup with Target cells in **sample A**.
5. Dispense 100 µl of **Target cell** suspension (**sample A**) to each assay well or tube.

Preparation of Effector cells

1. Prepare Effector cells at the appropriate concentration in **Wash Buffer**. For example, for a final Effector to **Target** ratio of 5:1, suspend Effector cells at 5×10^6 cells/ml.

Coincubation of Target and Effector cells

1. Add 100 µl of Effector cell suspension to each well (or tube) containing **Target cells** except at least two wells, and add 100 µl of Effector cell suspension to at least two wells which do not contain **Target cells**.
2. Add 100 µl of **Wash Buffer** to wells containing only **Targets** and only Effectors to bring all samples to a final volume of 200 µl.
3. Centrifuge all samples, carefully remove medium, and resuspend cell pellets in either 75µl **Substrate** or **Wash Buffer** for controls (absence of **Substrate** (for Tubes A and C below)). **For ADCC antibodies should be added at this time.**
4. Immediately after resuspension, pellet cells by brief centrifugation (Once the centrifuge reaches the speed normally used for cells, hold for 30 seconds and then stop the instrument.).
5. Incubate at 37°C for the desired time points. Since this assay detects dying cells rather than cells with irreversibly damaged plasma membranes, incubation times for a given cell system should be significantly shorter than with other

methodologies. For single point assays, the suggested coincubation time is 1 hour (enduser times range between 30 minutes and 2 hours). Substantially longer times are not recommended.

6. Wash each sample with 200 μ l **Wash Buffer**.
7. Resuspend each sample in **Wash Buffer**, transfer to flow cytometry tubes or leave in plates if a plate reader is to be used, and analyze by flow cytometry.

Summary of samples:

- A. Target cells**
 - A. Target cells** loaded with **TFL4** and **NFL1**
 - Ac. Target cells** loaded with **TFL4** only.
- B. Target cells (sample A) + Substrate**
- C. Effector cells**
- D. Effector cells + Substrate**
- E. Target cells(sample A) + Effector cells + Substrate** (multiple samples)

Flow Cytometry: the following channels are recommended: **PacBlue** for **NFL1**, **APC** for **TFL4**, and **FITC** for **substrate**. For instruments with other filter sets, please use settings consistent with the following excitation and emission peaks:


NFL1 λ_{ex} : 405 nm, λ_{em} : 450 nm

TFL4 λ_{ex} : 633 nm, λ_{em} : 657 nm

GS λ_{ex} : 488 nm, λ_{em} : 520 nm

For flow cytometers equipped with **Violet** (405), **Ar ion** (488 nm) and **red He-Ne** (633 nm) lasers:

1. Use samples **A** and **Ac** to eliminate target cells that have died prior to experiment. Using the **PacBlue** and **APC** channels, apply an **NFL1** gate to sample **A** as shown below and collect all cells contained therein.



2. Use samples **A** and **D** to initially set the **APC** and **FITC** channels, resp.: place the peak for cells from sample **A** near 10³ in the **APC** channel and the peak from sample **D** at ca. 10¹ in the **FITC** channel. Cells in sample B should then be at ca. 10³ in the **APC** and 10¹ in the **FITC** channels. (Note: Healthy (>95% viable by Trypan Blue) **TFL4**-labeled **Target cells** should appear as a single population. If not, first try decreasing labeling time to 15 minutes. If more than one population still exists, decrease **TFL4** concentration.)
3. Run remaining samples. 5,000-10,000 **Target cells** per sample should be collected for analysis.